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GÜVENLİLİK VE KLİNİK PERFORMANS ÖZETİ/SSCP **SUMMARY OF SAFETY AND CLINICAL PERFORMANCE/SSCP**

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This summary of safety and clinical performance (SSCP) aims to provide public access to an updated summary of key aspects of the safety and clinical performance of the device.

The SSCP is not intended to replace the instructions for use as the primary document to ensure the safe use of the device and is not intended to provide diagnostic or treatment recommendations to intended users or patients.

The following information is intended for users/healthcare professionals. This information is followed by a summary for patients.

1. CİHAZ TANIMI VE GENEL BİLGİLER

1. DEVICE DESCRIPTION AND GENERAL INFORMATION

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Strip; It is used especially as a suture line support material and tissue support - repair material in many surgical applications, including cardiovascular surgery, and is used to assist treatment by preventing tissue tearing at the suture line.

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1.6. Device Class: 2017/745/EU Annex-IX Conformity Assessment Based on Evaluation of Quality Management System and Technical Documentation, 2017/745/EU Annex-VIII Rule-8, Class III

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1.7. Year in which the First Certificate (CE) Covering the Device was Issued: The product received its first CE certificate in accordance with 93/42/EEC from the Approved Organization with code 2195 on 02.02.2014.

- 1.8. Varsa Yetkili Temsilci Adı ve SRN: --

1.8. Authorized Representative Name and SRN, if any: --

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1.9. Notified Body Name (notified body that will validate the SSCP) and Identification Number:--

2. CİHAZIN KULLANIM AMACI / 2. PURPOSE OF USE OF THE DEVICE

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2.1. Intended Use: PTFE (Polytetrafluoroethylene) Felt, Pledget, Strip; It is used especially as a suture line support material and tissue support - repair material in many surgical applications, including cardiovascular surgery, and is used to assist treatment by preventing tissue tearing at the suture line.

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2.2. Indication(s) and Target Population(s): Use as suture line support material and tissue support-repair material in many surgical applications, including cardiovascular surgery, and to prevent tissue tearing in the suture line. In cases where the surgical suture is likely to cause tearing in the tissue, there should be support under the suture. Various surgical suturing procedures such as vascular closure, septal repair, cardiac muscle closure, liver repair and valve suturing.

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2.3. Contraindications And/Or Limitations: Allergy to PTFE

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3.1. Description of the Device: PTFE (Polytetrafluoroethylene) Felt, Pledget, Strip; It is used especially as a suture line support material and tissue support - repair material in many surgical applications, including cardiovascular surgery, and is used to assist treatment by preventing tissue tearing at the suture line.

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Felt; It is used as tissue support material along with the suture line in tissue repair.

Strip; It is used as tissue support material along with the suture line in tissue repair.

Pledget; It is used as tissue support material along with the suture line in tissue repair.

PTFE (Polytetrafluoroethylene) Felt is used as a suture line support material in many surgical applications and prevents tissue tearing.

Felt, Strip and Pledget are made from the same material and are all included in the Felt product line. There is only a difference in size between them. Only the size of the products varies depending on the choice of the physician who will perform the operation and the intended use.

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3.2. A Reference to Previous Generation(s) or Variants, If Any, and an Explanation of Differences: Although our product has equivalent devices in the market, there are no previous or similar generations

Hazırlayan:	Onaylayan:
-------------	------------

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- Do not expose PTFE products to temperatures higher than 550 F (260 C). PTFE decomposes at high temperatures and forms highly toxic decomposition substances.
- Unused parts of the products should be thrown into the medical waste bin.

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4.3. Other Relevant Aspects of Safety, Including a Summary of Field Safety Corrective Actions (FSCA Including FSN), if applicable: Not available.

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5. SUMMARY OF CLINICAL EVALUATION AND POST-MARKET CLINICAL FOLLOW-UP (PMCF)

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5.1. Summary of Clinical Data Regarding the Equivalent Device, If Available:

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There are PTFE Teflon Felt Tissue Support Materials from many companies on the market. BARD Peripheral Vascular, Inc. Comparison with PTFE Teflon Felt of the company is given in the table below. The intended use, indications, contraindications and side effects of the products are the same. As seen below, the raw materials used, color and all other features are the same. Although there are minor differences, they are discussed in the clinical evaluation report.

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Warnings:

Do not expose PTFE products to temperatures higher than 550 F (260 C). PTFE decomposes at high temperatures and forms highly toxic decomposition pro-

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WARNINGS

1. This device is intended for single use only. If this medical device is to be re-used, medical devices (especially those with slits, joints and/or long and small lumens between components) are difficult or impossible to clean after an indefinite period of contact with tissues or body fluids with possible pyrogenic or microbial contamination. There is a risk of patient cross contamination. Residues of biological material may lead to contamination of the device with pyrogens or microorganisms that may lead to infectious complications.
2. Do not resterilize. The sterility of the product after resterilization is not guaranteed, as the level of potential pyrogenic or microbial contamination that could lead to infectious complications is uncertain. Cleaning, reprocessing, and/or resterilizing this medical device increases the possibility of the device not functioning properly due to components being adversely affected by thermal and/or mechanical changes.
3. Care should be taken with surgical tapes to ensure the integrity of the knots when used for permanent attachment.
4. Do not expose PTFE products to temperatures higher than 550°F (260°C). PTFE decomposes at high temperatures and forms highly toxic decomposition substances.

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ducts.

Unused parts of the products should be thrown into the medical waste bin.

Measures:

- The product sterilization bag must be opened in the operating room, which is a sterile environment, right before the operation and should not be used with bare hands in order to prevent contamination caused by the user.
- The product should only be used by surgeons specialized in relevant surgeries.
- The product size to be used in the operation may vary depending on the decision of the surgeon who will perform the operation.
- The surgeon should keep PTFE Felt in different sizes with him to use when a different size is required.
- Before use, it should be checked whether the package is damaged in order to ensure the sterility of the product.
- Contact of the product with any object should be avoided.
- The product should not be re-used under any circumstances.
- It should be transported and stored in environments without direct sunlight.
- The user manual should be read carefully before using the product.
- Torn and damaged packaged products are not sterile and should not be used.
- Transportation and storage should be done in a way that the packages will not be dama-

MEASURES

1. If the packaging is not open or damaged, it is STERILE. It is for single use only.
2. These products are sterilized with ethylene-oxide. DO NOT RE-STERILIZE.
3. In order not to damage this braided product, care must be taken when pressing the adhesive tapes.
4. When surgical adhesive tapes are used as placement material, it is recommended to use tapered and non-cutting needles. Thus, the fibers will not be cut during sewing.
5. Avoid touching the felt surface too much before use.
6. When cutting felt into smaller pieces it is important to use a sharp cutting device or scissors with sharp blades to ensure that the felt fibers are cut smoothly and are not pulled apart.
7. After use, this product may contain potential biohazards. Use and dispose of the product only in accordance with established medical practice and in accordance with applicable laws and regulations.

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	ged. ➤ Heavy materials should not be placed on the boxes. ➤ Expired products should not be used. ➤ This product should not be re-sterilized under any circumstances. ➤ The entire product or the remainder of the product should not be re-sterilized by the purchaser. ➤ Despite this, the manufacturer and supplier do not accept any responsibility when the product is sterilized again. ➤ If there is a foreign substance in the product and it is used, it may cause infection and toxic effects for the patient. If there is foreign material in the product, the product should not be used. ➤ If the product deteriorates during use and this product is used, the expected performance cannot be obtained from the product. It causes the treatment to be prolonged. Spoiled products should not be used. ➤ If force is applied during product use and the product deteriorates and is used, the product cannot be expected to perform as expected. The product should not be subjected to force or damaged. ➤ When cutting felt into smaller pieces, it is important to use a sharp cutting device (scalpel) or scissors with sharp blades to ensure that the felt fibers are cut smoothly and are not pulled apart.	
Kontraendikasyonlar Contraindications	PTFE karşı alerji Allergy to PTFE	PTFE karşı alerji Allergy to PTFE

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5.2. Summary of Clinical Data Obtained from Studies Conducted Before CE Marking of the Device, If Any: --

5.3. Summary of Clinical Data Obtained from Other Sources, If Any: In line with data obtained from after-sales surveillance studies, risk management studies, literature and other databases, it has been

Hazırlayan:	Onaylayan:
-------------	------------

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determined that the clinical benefit of our product is much greater than its risk.

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5.4. An Overview of Clinical Performance and Safety:

PTFE (Polytetrafluoroethylene) Felt, Pledget, Strip; It is used as a suture line support material and tissue support-repair material in many surgical applications, including cardiovascular surgery, and is used to assist treatment by preventing tissue tearing at the suture line. Indications; It is used as a suture line support material and tissue support-repair material in many surgical applications, including cardiovascular surgery, and to prevent tissue tearing at the suture line. In cases where the surgical suture is likely to cause tearing in the tissue, there should be support under the suture. Various surgical suturing procedures such as vascular closure, septal repair, cardiac muscle closure, liver repair and valve suturing.

A system has been established and monitored to track recalls and adverse events of equivalent/similar devices.

This clinical evaluation confirmed compliance with the General Safety and Performance Requirements in MDR 2017/745.

All data support that PTFE Felt, Pledget, Strip products perform acceptably and as expected, with no safety concerns noted and the overall residual risk being acceptable when weighed against the benefits to the patient.

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5.5. Ongoing or Planned Post-Market Clinical Follow-up: In accordance with MDR Annex XIV Part B, CE-marked and In accordance with the MDCG 2020-7 and MDCG 2020-8 Guidelines, for the purpose of proactively collecting clinical data obtained from the use on or within the human body of a device placed on the market or put into service within its intended use, as specified in the relevant conformity assessment procedure. PMCF work is carried out in line with the Protocol and Plan created in accordance with the plan.

6. OLASI TANI VEYA TEDAVİ ALTERNATİFLERİ / POSSIBLE DIAGNOSIS OR TREATMENT ALTERNATIVES

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Our product is used to support the suture line. Our product consists of PTFE. Publications using Dacron as an alternative raw material can be accessed through literature review. There are also publications comparing alternative adhesives (bioglue) with PTFE Felt-Strip-Pledget. In general, it can be seen that adhesives cause different risks.

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The device can only be applied by a doctor who is fully responsible for the surgical and post-operative procedures of the system and can manage any possible problems. Correct selection of PTFE Felt application area, thickness and size for each patient is extremely important for the successful implementation of the procedure. Therefore, the device must be used by a healthcare professional / doctor (operating surgeon).

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Sıra No Sequence No.	Doküman No Document No	Doküman Adı Document Name
1)	EN ISO 13485:2016/A11	Medical devices - Quality management systems
2)	EN ISO 14971:2019/A11	Medical devices - Application of risk management to medical devices
3)	EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
4)	EN ISO 10993-3	Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity
5)	EN ISO 10993-5	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
6)	EN ISO 10993-4	Biological evaluation of medical devices - Part 4: Selection of tests for interactions with blood
7)	EN ISO 10993-6	Biological evaluation of medical devices - Part 6: Tests for local effects after implantation

Hazırlayan:	Onaylayan:
-------------	------------

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8)	EN ISO 10993-7+2008/A1	Biological evaluation of medical devices -- Part 7: Ethylene oxide sterilization residuals
9)	EN ISO 10993-10	Biological evaluation of medical devices -- Part 10: Tests for irritation and skin sensitization
10)	EN ISO 10993-11	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
11)	EN ISO 10993-12	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
12)	EN ISO 10993-17	Biological evaluation of medical devices — Part 17: Establishment of allowable limits for leachable substances
13)	EN ISO 10993-18	Biological evaluation of medical devices - Part 18: Chemical characterization of medical device materials within a risk management process
14)	EN 15223-1	Tıbbi Cihazların Etiketlenmesinde Kullanılan Semboller
15)	EN 62366-1:2015/A1	Medical devices - Part 1: Application of usability engineering to medical devices
16)	EN ISO 20417	Medical devices - Information to be supplied by the manufacturer
17)	EN ISO 11607-1	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
18)	EN ISO 11607-2	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes
19)	EN 865-5	Packaging for terminally sterilized medical devices - Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods
20)	EN ISO 14937	Sterilization of health care products - General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices
21)	EN ISO 11138-2	Sterilization of health care products - Biological indicators - Part 2: Biological indicators for ethylene oxide sterilization processes
22)	EN ISO 11737-1	Sterilization of medical devices - Microbiological methods - Part 1: Determination of a population of microorganisms on products
23)	EN ISO 11737-2	Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
24)	EN ISO 14644-1	Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness by particle concentration
25)	EN ISO 14644-2	Cleanrooms and associated controlled environments - Part 2: Monitoring to provide evidence of cleanroom performance related to air cleanliness by particle concentration
26)	EN ISO 14644-3	Cleanrooms and associated controlled environments - Part 3: Test methods
27)	EN ISO 14644-4	Cleanrooms and associated controlled environments - Part 4: Design, construction and start-up
28)	EN ISO 14644-5	Cleanrooms and associated controlled environments - Part 5: Operations
29)	EN ISO 14644-7	Cleanrooms and associated controlled environments — Part 7: Separative devices (clean air hoods, gloveboxes, isolators and mini-environments)
30)	EN ISO 14644-8	Cleanrooms and associated controlled environments — Part 8: Classification of air cleanliness by chemical concentration (ACC)
31)	EN 17141	Cleanrooms and associated controlled environments - Biocontamination control
32)	EN ISO 11140-1	Sterilization of health care products - Chemical indicators - Part 1: General requirements
33)	TS EN 556-1	Sterilisation of medical devices - Requirements for medical devices to be designated "Sterile" - Part 1: Requirements for terminally sterilised medical devices
34)	EN ISO 11135:2014/A1	Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices
35)	EN ISO 14630	Non-active surgical implants - General requirements
36)	ASTM F 1980	Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices
37)	ASTM D5729	Standard Test Method for Thickness of Nonwoven Fabrics
38)	ASTM D3776	Standard Test Methods for Mass Per Unit Area (Weight) of Fabric
39)	ASTM D5034	Standard Test Method for Breaking Strength and Elongation of Textile Fabrics (Grab Test)

Hazırlayan:	Onaylayan:
-------------	------------

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40)	ASTM D3786	Standard Test Method for Bursting Strength of Textile Fabrics—Diaphragm Bursting Strength Tester Method
-----	------------	---

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A summary of the safety and clinical performance of the device for patients is provided below.

GÜVENLİLİK VE KLİNİK PERFORMANS ÖZETİ/SSCP
SUMMARY OF SAFETY AND CLINICAL PERFORMANCE/SSCP

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This summary of safety and clinical performance (SSCP) aims to provide public access to an updated summary of key aspects of the safety and clinical performance of the device. The information presented below is intended for patients or laypeople. A more comprehensive summary of the safety and clinical performance of the device for healthcare professionals is included in the first section of this document.

Hazırlayan:	Onaylayan:
-------------	------------

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SSCP is not intended to give general advice about treating a medical condition. [If you have questions about your medical condition or the use of the device in your condition, please contact your healthcare professional.](#) This SSCP is [not intended to replace an implant card or instructions for use to provide information on the safe use of the device.](#)

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1. DEVICE IDENTIFICATION AND GENERAL INFORMATION

- Trade Name of the Device: PTFE TEFLON FELT-PLEDGET-STRIP
- Manufacturer Name and Address: GEMED Sağlık Ürünleri Sanayi Ve Ticaret Ltd. Ltd. - Hürriyet, Adnan Kahveci Cd. No:27/B, 34876 Kartal/ İstanbul/ Türkiye
- Basic UDI-DI: 869907384TD01FQH
- Year in which the CE Mark was first given to the device: The product received its first CE certificate in accordance with 93/42/EEC from the Approved Body coded 2195 on 02.02.2014.

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1. INTENDED USE OF THE DEVICE

- Intended Use: PTFE (Polytetrafluoroethylene) Felt, Pledget, Strip; is used in many surgical applications, including cardiovascular surgery, especially as a suture line support material and tissue support - repair material, and is used to help treatment by preventing tissue tearing in the suture line.
- Felt; is used as a tissue support material together with the suture line in tissue repair.
- Strip; is used as a tissue support material together with the suture line in tissue repair.
- Pledget; is used as a tissue support material together with the suture line in tissue repair.
- PTFE (Polytetrafluoroethylene) Felt is used in many surgical applications especially as a suture line support material and prevents tissue tearing.
- Felt, Strip and Pledget are made of the same material and are all included in the Felt product group. There is only a difference in size between them. Only the size of the products changes according to the option of the physician who will perform the operation and the purpose of use.
- Indications and Target Patient Groups: In many surgical applications including cardiovascular surgery, it is used especially as suture line support material and tissue support - repair material and prevents tissue tearing in the suture line. In cases where there is a possibility of surgical suture causing tissue tearing, it provides support under the suture site. Various surgical suture procedures such as vascular closure, septal repair, heart muscle closure, liver repair and valve suture.
- It is used in patients who need to use it as a support material in surgical interventions and especially in suture line operations to provide non-absorbent tissue protection during surgical suture procedures. These are patients who require surgical interventions including cardiac surgery and who require support material for the suture line. It can be used in women or men according to the intended use.
- Contraindications: Allergy to PTFE

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3. DEVICE DESCRIPTION

- Device Description and Materials/Substances in Contact with Patient Tissues: PTFE (Polytetrafluoroethylene) Felt, Pledget, Strip; It is used in many surgical applications including cardiovascular surgery, especially as suture line support material and tissue support - repair material and is used to help treatment by preventing tissue tearing in the suture line.
- Felt; It is used as tissue support material together with suture line in tissue repair.
- Strip; It is used as tissue support material together with suture line in tissue repair.
- Pledget; It is used as tissue support material together with suture line in tissue repair.
- PTFE (Polytetrafluoroethylene) Felt is used in many surgical applications especially as suture line support material and prevents tissue tearing.
- Felt, Strip and Pledget are made of the same material and are all included in the Felt product group. There is only a difference in size between them. Only the size of the products changes according to the choice of the physician who will perform the operation and the purpose of use.
- The entire product comes into contact with the heart and remains in the patient's body for the patient's life. In this direction, there are general explanations of biological evaluation.
- Information about Medical Products in the Device, if any: --
- Explanation of How the Device Achieves Its Intended Effect Mechanism: PTFE (Polytetrafluoroethylene) Felt is used as a suture line support material in many surgical applications and prevents tissue tearing. In cases where there is a possibility that the surgical suture will cause tissue tearing, it provides support under the suture. While the product shows the specified effect, it does not show any chemical or biological reaction results. Therefore, it physically acts in accordance with the purpose of use with the implantation.
- Explanation of Accessories, if any: The device is offered to the market in single-use, sterile packaging and does not have any accessories.

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4. RISKS AND WARNINGS

If you have any concerns about side effects or risks related to the device or its use, contact your healthcare professional. This document is not intended to replace consulting your healthcare professional when necessary.

o How Potential Risks Are Controlled or Managed:

o Residual Risks and Unwanted Effects: The product may have an allergic effect when applied to patients with PTFE allergy. Infections that may occur in any surgical operation may occur.

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- The product sterilization bag must be opened immediately before the operation in a sterile operating room and should not be used with bare hands to prevent contamination from the user.
- The product should only be used by surgeons who are experts in the relevant operations.
- The product size to be used in the operation may vary depending on the decision of the surgeon who will perform the operation.
- The surgeon should have PTFE Felt in different sizes with him/her in case a different size is required.
- Before use, the product should be checked for damage to its packaging in terms of sterility.
- The product should not come into contact with any object.
- The product should not be reused under any circumstances.
- It should be transported and stored in environments that are not exposed to direct sunlight.
- The user manual should be read carefully before using the product.
- Products with torn and damaged packages are not sterile and should not be used.
- Transportation and storage should be done in a way that the packages will not be damaged.
- Heavy materials should not be placed on the boxes.
- Products that have expired should not be used.
- This product should not be re-sterilized under any circumstances.
- The entire product or the remaining product should not be re-sterilized by the purchaser.
- Despite this, the manufacturer and supplier assume no responsibility when the product is re-sterilized.
- If there is a foreign substance in the product during use and it is used, it will cause infection and toxic effects for the patient. The product should not be used if there is a foreign substance in the product.
- If the product is deteriorated during use and this product is used, the expected performance cannot be obtained from the product. It causes the treatment to be prolonged. Deteriorated products should not be used.
- If force is applied during product use and the product deteriorates and is used, the product cannot be expected to perform as expected. Force should not be applied to the product, it should not deteriorate.
- When cutting felt into smaller pieces, it is important to use a sharp cutting tool (scalpel) or scissors with sharp blades to ensure that the felt fibers are cut properly and not pulled apart.

Hazırlayan:	Onaylayan:
-------------	------------

**GÜVENLİLİK VE KLİNİK PERFORMANS
ÖZETİ RAPORU - SSCP
SAFETY AND CLINICAL PERFORMANCE
SUMMARY REPORT - SSCP**

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• Do not expose PTFE products to temperatures greater than 550 degrees F (260 degrees C). PTFE decomposes at high temperatures and forms highly toxic decomposition products.

• Unused product should be disposed of in a medical waste container.

o Summary of Field Safety Corrective Action (FSCA Including FSN), if applicable: Not available.

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5. CLINICAL EVALUATION SUMMARY AND CLINICAL FOLLOW-UP AFTER MARKETING

Clinical History of the Device:

Our product has been on the market for many years. Again, products with the same features and raw materials have been used in operations for many years by many manufacturers other than us. PTFE (Polytetrafluoroethylene) Felt, Pledget, Strip; It is used in many surgical applications including cardiovascular surgery, especially as suture line support material and tissue support - repair material and is used to help treatment by preventing tissue tearing in the suture line. Indications; It is used in many surgical applications including cardiovascular surgery, especially as suture line support material and tissue support - repair material and to prevent tissue tearing in the suture line. In cases where there is a possibility that the surgical suture will cause tissue tearing, support under the suture site. Various surgical suture procedures such as vascular closure, septal repair, heart muscle closure, liver repair and valve suture. A system has been established and is being monitored to track recalls and adverse events of equivalent/similar devices.

This clinical evaluation has confirmed compliance with the General Safety and Performance Requirements given in MDR 2017/745.

All data support that PTFE Felt, Pledget, Strip products perform acceptably and as expected, no safety issues have been noted and the overall residual risk is acceptable when weighed against the benefits to the patient.

Clinical Evidence for CE Marking:

PTFE (Polytetrafluoroethylene) Felt, Pledget, Strip; is used in many surgical applications including cardiovascular surgery, especially as a suture line support material and tissue support - repair material and is used to prevent tissue tearing at the suture line and to assist treatment. Indications; It is used in many surgical applications including cardiovascular surgery, especially as a suture line support material and tissue support - repair material and to prevent tissue tearing at the suture line. Support under the suture site in cases where there is a possibility of tissue tearing due to surgical suture. Various surgical suture procedures such as vascular closure, septal repair, cardiac muscle closure, liver repair and valve suture.

Hazırlayan:	Onaylayan:
-------------	------------

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A system has been established and is being monitored to track recalls and adverse events related to equivalent/similar devices.

This clinical evaluation has confirmed compliance with the General Safety and Performance Requirements given in MDR 2017/745.

All data support that PTFE Felt, Pledget, Strip products perform acceptably and as expected, no safety issues have been noted and the overall residual risk is acceptable when weighed against the benefits provided to the patient.

Safety:

Our product has been on the market for many years. Again, products with the same properties and raw materials have been used in operations by many manufacturers other than us for many years. PTFE (Polytetrafluoroethylene) Felt, Pledget, Strip; It is used in many surgical applications, including cardiovascular surgery, especially as a suture line support material and tissue support - repair material, and is used to help treatment by preventing tissue tearing in the suture line. Indications; It is used in many surgical applications, including cardiovascular surgery, especially as a suture line support material and tissue support - repair material, and to prevent tissue tearing in the suture line. In cases where there is a possibility that the surgical suture will cause tissue tearing, support under the suture site. Various surgical suture procedures such as vascular closure, septal repair, heart muscle closure, liver repair and valve suture.

A system has been established and is being monitored to track recalls and adverse events related to equivalent/similar devices.

This clinical evaluation has confirmed compliance with the General Safety and Performance Requirements given in MDR 2017/745.

o All data support that PTFE Felt, Pledget, Strip products perform acceptably and as expected with no safety issues noted and the overall residual risk is acceptable when weighed against the benefits to the patient.

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6. POSSIBLE DIAGNOSIS OR TREATMENT ALTERNATIVES

When considering alternative treatments, it is recommended that you consult your healthcare professional who can take your individual situation into consideration.

o General Description of Treatment Alternatives: Our product is used to support the suture line. Our product consists of PTFE. Publications using Dacron as an alternative raw material can be found by lite-

Hazırlayan:	Onaylayan:
-------------	------------

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ature search. There are also publications comparing alternative adhesives (bioglue) with PTFE Felt-Strip-Pledget. In general, it can be seen that adhesives cause different risks.

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The products are used by healthcare professionals. It is not possible for the patient to use the product.

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